

# SASS<sup>®</sup> 3010

## Particle Extractor



Filter-to-liquid particle extraction, simplified.

### Features

- Fast and efficient extraction from electret filters
- Extraction efficiency range 70 – 80%, typ.

### Application Areas

- Pharmaceutical
- Medical facilities
- Public health
- Clean rooms
- Military
- Food processing
- UAVs
- Agriculture
- Indoor air quality
- Environmental
- Homeland security

**The SASS<sup>®</sup> 3010** is used to extract and transfer aerosols from Research International's electret filters to a small fluid volume for analysis.

**Ease of use** Captured particulates can be difficult to remove because induced dipole fields create a strong holding force and must be neutralized. Once particulates have been released, they must then be removed from within the fibrous filter matrix and collected in a small amount of sample fluid. These processes are efficiently performed in a matter of 1–2 minutes using the SASS 3010 Particle Extractor.

**Tested for efficiency** To test extraction efficiency, several electret filters were used to collect



airborne fluorescent polystyrene beads of 1.8 microns diameter. Each filter was operated for a period of 10 minutes. After the collection phase was completed, the filters were mounted in the SASS 3010 and captured beads transferred to 5 ml of extraction buffer. Extraction efficiencies were then determined using fluorometric assay methods.

**Results** It was found that an average recovery of 77% was achieved. A second extraction with an additional 5 ml of extraction fluid resulted in recovery of another 17% of the embedded beads, while two more 5 ml extractions resulted in small 4.5% and 1.5% additions to the total number of beads recovered, respectively.

## SASS® 3010 Specifications

<b>Filter compatibility</b>	For use with SASS 3100 and SASS 4100 filters.
<b>Extraction method</b>	Acoustic vibration of the fluid-saturated filter is followed by counter-flow discharge of the suspended aerosol particles.
<b>Extraction efficiency</b>	70-80% typical
<b>Carry-over</b>	1.1% with dry wiping, and 0.01% to 0.1% with a 5 ml flush. Additional flushes will reduce carry-over further.
<b>Extraction time</b>	1 to 2 minutes, typical, with a flush cycle.
<b>Sample fluid storage</b>	The extraction fluid bottle is also used for fluid sample storage upon completion of extraction.
<b>Physical size</b>	<ul style="list-style-type: none"> <li>• Body: 10.2 cm (W) x 13.4 cm (D) x 14.5 cm (H)</li> <li>• A 7.8 cm-high plunger protrudes from the extractor's top surface.</li> </ul>
<b>Weight</b>	800 grams
<b>Electrical power</b>	Two size "D" primary batteries. Optional power supply available, purchased separately.
<b>Operating temperature</b>	0°C to 70°C
<b>Extraction fluid</b>	A pre-filled dropper bottle provides enough buffered extraction fluid to make a 5 ml sample. Other fill levels to 10ml available on user request. Purchased separately.



*Research International reserves the right to change specifications without prior notice.*



### U.S. Headquarters Office

Research International, Inc.  
17161 Beaton Rd. S.E.  
Monroe, WA 98272-1034 USA

Tel: 1.800.927.7831  
info@resrchintl.com  
www.resrchintl.com

### U.S. East Coast Office

Mr. Jon Tobelmann  
1.703.625.8381  
jontobelmann@resrchintl.com

### International Offices

Please contact the U.S. Headquarters office to locate a representative in your region.